



Your Veterinary Microbiology Specialists



Welcome to Hardy Diagnostics!

Why is Hardy Diagnostics the right choice for your lab?

* There is no "middle man". You will receive your supplies directly from Hardy Diagnostics ISO certified factory. This eliminates the problems associated with a "middle man" and ensures a long product shelf life.



* We're available. Feel free to contact us anytime to discuss your technical needs and questions. Our focus on customer service is noted as being the best in the industry!

Technical Support:

800-266-2222, Extension 5658

TechService@HardyDiagnostics.com

Customer Service:

800-266-2222

CustService@HardyDiagnostics.com

- * We offer a huge selection. You can choose from a complete selection of microbiology and laboratory products—over 13,000 of them!-including many "hard to find" items.
- * **Speedy service.** Thanks to Hardy Diagnostics' network of seven distribution centers across the U.S. (California, Washington, Utah, Arizona, Texas, Ohio, Florida), your shipments are delivered quickly, usually the next day.

We are standing by to fulfill your orders and we are ready to assist you with any questions that you might have.

Sincerely,

Jay Hardy, MT(ASCP), SM(NRCM)

President



ISO 13485 FM 572526

Table Of Contents

General Lab Supplies

el la Neel Ive I	
Gloves: Latex, Nitrile and Vinyl	
Gold Gloves1	
Precision® Gloves1	
VPF Gloves1	
Autoclave/Biohazard Bags2)
Autoclave/Biohazard Bags2	
Autoclave Bags2	
FreshClave [™]	
Biological Sterility Indicator2	
Biohazard Bag3	
Sharps Container3)
Strar ps Cortainer)
Laboratory Attire4	ļ
ComfortPRO™ Lab Coats4	ļ
Face Shield4	ļ
Lab Apron4	ļ
Isolation Mask with Ear Loops4	1
Bouffant Cap4	ļ
D: . / T . T .	
Pipets/Test Tubes	
Culture Tube5	
Centrifuge Tube5	
Capillary Tube5)
Microscope Supplies6	5
ProSlide [™] , Microscope Slides6	
Cover Glass6	
Immersion Oils6	ó
Lens Paper6	
Microscope Lens Cleaner6	
Bibulous Paper	
Slide Mailer, Cardboard	
Slide Mailer, Polypropylene7	7
Slide Storage Box	
S .	
Swabs & Applicators8	
Cotton Swabs, Sterile8	
Transporter Swabs8	
Tongue Blade/Sampling Stick8	
Dacron® Swabs, Sterile8	3
Specimen Containers)
SureTight™ Specimen Sample Cups	
Whirl-Pak® Write-On Specimen Bags	
·	
Disinfectants/Cleansers, Sterilants10	
Alpet D210)
CitriFoam® Disinfectant/Cleaner10	
CitriGuard® Disinfectant/Cleaner10)

DisCide® Towelettes	10
Vavicide®	
Hand Soaps, Lotions, Sanitizers	11
AloeGuard® Antimicrobial Soap	11
AloeSoothe™ Moisturizing Lotion	11
AloeSafe [™] Hand Sanitizer Gel	11
Smart-San [™] Healthcare Personnel Handwash	11

Microbiology

Bovine Mastitis	12
Bovine Blood Agar	
Bovine Blood Agar Plate	12
Bovine Blood Agar Plate with Esculin	12
Bovine Blood Agar with Esculin, Modified	
Bovine Selective Strep Agar	12
Bovine Selective Staph Agar	
Mastitis Biplate	12
Mastitis Triplate	12
Mastitis Quadplate	
InPouch TF-Bovine	12
Urinalysis	13
Urine Dipstick Test	
Castile Soap Towelettes	
UriStain [™]	
SureTight [™] Sample Cups	
Refractometer	13
Urine Culture System	
HUrBi [™]	
Blood Agar, MacConkey Agar Biplate	
Dip-Slides	15
Skin Culture Systems (Dermatophytes)	
Sab-Duet [™] Biplate	
DTM (Dermatophyte Test Medium)	
Derm-Duet [™]	14
HardCHROM [™] Products	
HardyCHROM [™] Staph aureus	
HardyCHROM™ MRSA	
HardyCHROM [™] Candida	15
Culture Media	
Bile Esculin Agar	
Blood Agar, 5%/MacConkey	
Blood Agar, 5% Sheep in TSA Base	16

Brucella Agar with Hemin		McFarland Standard	22
and Vitamin K /LKV	16	Wickerham Card	
Campy CVA Agar		Sterile Swabs	
Chocolate Agar		Sterne Swabs	
Columbia Blood Agar/MacConkey		Loops and Innoculation	23
EnteroScreen4 [™]		Disposable SpeedStreaks [™]	
GN Broth		Reusable Calibrated Nichrome Loop	
HE Agar		Bacti-Cinerator™ III	
LIA (Lysine Iron Agar) Slant		Spreaders	
MacConkey Agar		5p1eade13	
Mueller Hinton Agar		Staining Supplies	24
Mueller Hinton Agar w/Bld		Gram Stain Advanced™ Kit	24
Mueller Hinton Agar		Erada-Stain Stain Remover	
Mycoplasma Agar		Coplin Jar	
Rose Agar (CNA and PEA combination)		Slide Staining Rack	
		Q-Slide Gram (Staining Control)	ر∠
Thioglycollate with H and K		Stat Stain, Rapid Wright-Giemsa	25
SabDex C&G/DTM Biplate		Fecal Fat StainAcid Fast Bacterium Stain	
SabDex (Sabouraud Dextrose) Agar Slant			
Selenite Cystine Broth		Mayer's Albumin Fecal Occult Blood: QuickVue® iFOB	
Thio w/o Indicator		recai Occuit Biood: Quickvue* irOB	20
TSI Agar Slant			
Urea Agar Slant	1 /	Parasitology	
Reagents	18	Diagnostics Tests	26
Bacitracin HardyDisk [™]	18	Para-TECT [™] Giardia Antigen Detection	
Catalase Reagent		InPouch™ TF-Bovine	
CoaguStaph [™]		InPouch™ TF-Feline	
OxiDrops [™]		PARA-TECT™ Cryptosporidium Slide	
India Ink, Dropper Bottle		ImmunoCard Stat! Crypto/Giardia	
Indole		mmunocura stat. Cryptor diaraia	
Nitrate Reagent A Sulfanilic Acid	18	Control Slides	27
Nitrate Reagent B Sulfanilic Acid		Cryptosporidium Slide	
Optochin, HardyDisk™		Intestinal Protozoa	27
OxiStrips [™]		Trestriar i Totozoa	
Potassium Hydroxide		Concentration	
PYR Test Kit		ParaCon [™] Fecal Concentration Kit	28
Voges-Proskauer Reagent A	18	raiacon recarconcentration (it	20
Voges-Proskauer Reagent B		Transport	
3		Stool Transport Vial	28
RapidTests	19	Stool Hallsport vial	20
StaphTEX™ Blue Kit		Stains	
StrepPro [™] Grouping Kit		Trichrome Stain	28
MRSA Latex Test for PBP2		Trichrome Stain Starter Kit	
		10% Formalin	
Quality Control Microorganisms	19	Modified (Cu) PVA Fixative	
		(/	
Susceptibility Testing		Equipment	
Antibiotic Disks: HardyDisk™ AST		Equipment	
Mueller Hinton Agar Plates		Exact-Temp [™] Thermometers	29
Saline, 0.85%		Incubator	
Disk Dispensers	22		



HUGO[™] is free for all Hardy Diagnostics Customers!





HUGO[™] version 7 is a unique software program from Hardy Diagnostics that serves as an encyclopedia and technical manual for microbiologists.

This software is offered free to microbiologists as a valuable reference source. It will provide quick and easy access to important up-to-date information regarding culture media, procedures, quality control, safety, antimicrobials, agency regulations, new products, rapid methods, and other topics of interest to microbiologists.

General Lab Supplies Gloves: Latex, Nitrile and Vinyl



Gold Gloves

Powder-free, cream colored, latex gloves.

- ✔ Powder-free
- ✓ Eliminates powder related contamination in critical exam or in clean room environments.
- ✓ Full textured surface provides superior grip in wet or dry applications.
- ✓ Smooth interior finish for easy donning and comfortable wear.

10 boxes/case

Extra Small	GLD260
Small	GLD262
Medium	GLD265
Large	GLD266
Extra Large	GLD268



Precision® Gloves

Powder-free, violet colored nitrile, latex-free, gloves.

- Ultra soft and super stretchy for a very comfortable fit and to reduce hand fatigue.
- ✓ 100% latex-free eliminates potential allergic reactions to latex protein.
- Textured surface provides strong grip for both wet and dry applications.
- Superior puncture, tear, and chemical resistance with excellent dexterity.

10 boxes/case

Extra Small	PCS770
Small	PCS772
Medium	PCS775
Large	PCS776
Extra Large	



VPF (Vinyl) Gloves

Powder-free, translucent, vinyl gloves.

- √ 100% latex free and powder-free minimizes the risk of latex protein and chemical residue-related allergies.
- ✓ Eliminate powder-related contamination in critical exam or in clean room environments.
- ✓ Very soft and flexible to fit snugly and comfortably.
- ✓ Strong and durable for superior protection.

10 boxes/case

Extra Small	VPF230
Small	VPF232
Medium	VPF235
Large	VPF236
Extra Large	VPF238

General Lab Supplies Autoclave Bags



Autoclave Bags

Made of tough, thick polypropylene, doubleseamed bottom, solid red bag with black biohazard warnings

in English and Spanish. Bag is rated to withstand heat up to 300°F. Heat indicator strip turns black when exposed to steam, OSHA approved labeling.

8x12 inches, 2mil thickness 400/case	AB5819
12x24 inches, 2mil thickness 500/case	AB5010
14x19 inches, 2mil thickness 500/case	AB5012
19x24 inches, 2mil thickness 200/case	AB5013
25x30 inches, 1.75 mil thickness 200/case	AB5026
25x35 inches, 2mil thickness	AB5237

FreshClave[™] Autoclave Deodorant

Simple and convenient method of eliminating the unpleasant odors associated with disposal of laboratory waste. Presented in a soft gelatin capsule containing a concentrated liquid fragrance that melts in the heat of autoclave cycles: provides a comfortable, fragrant work environment. One bead will deodorize up to 10 cubic feet. 100 beads/bottle.



Apple scented	FC604
Lemon scented	FC605
Peach scented	FC603



SporAmpule® Liquid, Cat. no. BST20

Biological Indicator

An error-free, BI (Biological Indicator) system certified for use with steam sterilization processes. After autoclaving with a normal load, the SporAmpule® is incubated at 56°C for 48 hours. A purple to yellow color change and visible turbidity indicate sterilization failure, by SPS Medical.

20 per package	BST20
100 per package	BST100

General Lab Supplies Biohazard Bags



Biohazard Bag

For collection of biohazard waste, not for autoclaving. Made of polyethylene, solid red with black biohazard warnings in English and Spanish. Labeled with international biohazard symbol and meets OSHA requirements: passed the 165gm dart impact test.

Red, 11x14 inches, 1.2mil thickness, Polyethylene,1000/boxBB	2214
Red, 23x23 inches,1.5mil thickness, Polyethylene 500/boxB	B116
Red, 23x23 inches, 1.2mil thickness, Polyethylene, 500/boxBB	2021
Red, 24x32 inches, 1.2mil thickness, Polyethylene, 250/boxBB	2192
Red, 25x35 inches, 2.25mil thickness, Polyethylene, 200/boxBB	2305



General Lab Supplies Laboratory Attire



ComfortPRO™ LabCoats

For safety, style, and comfort! These affordable, disposable lab coats are made of tri-laminated, spunbound material (SMS) which is designed to be fully breathable and splash resistant. Our lab coats also feature breast and hip pockets, snap front buttons, knit cuffs and collars for added comfort.

Size	White	Blue
S	36103WS	36103BS
M	36103WM	36103BM
L	36103WL	36103BL
XL	36103WXL	36103BXL
XXI	36103WXXI	36103BXXI

Face Shield



Complete head band and window. Comfortable plastic headband that adjusts on the side and across the top, provides full face protection as required by OSHA, window is thick durable acetate, the

window can be flipped up when not in use.

Face Shield, with flip-up window, complete33007060

Face Shield replacement window only33007071



Lab Apron Personal, protective, disposable apron made of polyethylene, 2-1/4 mil. thick. Perfect for processing specimens where splashing may occur.

Plastic, fluid proof, 43x59 inches, 20 per case......8571



Isolation Mask with **Ear Loops**

Latex free with convenient earloops sonically welded to the outside of the mask. Soft, cool, and

breathable while meeting the ASTM F200-01 testing requirements of submicron filtration above 98%.

500 per case......15100




Bouffant Cap

Consistency in sizing, strength, and comfort. Made of polypropylene. Elastic band is sized for a secure. comfortable fit.


500 per case......MW2111

General Lab Supplies Pipets/Test Tubes



Features

- Non-toxic and chemically inert to biological Produced from a custom grade of lowfluids and most acids
- Can be heat sealed to use for liquid storage or transport, temperature range: -196°C to 90°C
- Pipets can be frozen in liquid nitrogen



Culture Tube

Disposable borosilicate glass culture. These tubes are chemically clean and well suited for general use in blood banks, hematology, chemistry

and all areas of the laboratory. Manufactured in compliance with strict specifications, with consistent size and weight and properly rounded bottoms.

Borosilicate glass, round bottom:

10x75mm,1000/box	735001075
12x75mm,1000/box	735001275
13x100mm,1000/box	7350013100
15x85mm,1000/box	735001585
16x100mm,1000/box	7350016100
16x125mm,1000/box	7350016125
16x150mm,1000/box	7350016150

- density polyethylene (LDPE) to make them non-slippery
- Easy-grip with gloves
- Custom packaging
- · Can be gas (EtO) sterilized

Centrifuge

Tube

Polypropylene, molded graduations, with screw cap. 15ml, by Globe Scientific.




500 per case6	275
250 per bag, 1000 per case6	261

Capillary Tube

Heparinized for hematocrits, contains heparin, red tip, 1.1mm inside diameter. Red tip, 75mm long, by Chase Instruments. 200 per vial......2501

General Lab Supplies Microscope Supplies



ProSlide[™], Microscope Slides

Guaranteed non-sticking, ideally suited for hematology differential smears, pathology tissue mounting and all general microscopy. Slides are manufactured from highest purity, corrosion-resistant clear sheet glass fully annealed to minimize chipping. Edges are precision ground for smoothness.

3x1 inch x 1mm, 72/box.	
Frosted, economy	PF72E
Frosted, premium	PF72P
Plain, economy	PP72E
Plain, premium	PP72P

Cover Glass

High quality, borosilicate cover glass is non-corrosive and non-fogging. Devoid of scratches, bubbles and striations, its flatness provides a uniform, distortion-free surface. A desiccant is included so that moisture doesn't degrade the cover glass. Convenient plastic storage container allows for easy stacking on the shelf. 1oz box.

18x18mm	1818HD
22x22mm	2222HD
22x40mm	2240HD
24x40mm	2440HD
24x50mm	2450HD



Immersion Oils

Type A, low viscosity	
4oz	Z85
20ml (with glass applicator)	Z95
Type B, high viscosity	
4oz	Z86
20ml (with glass applicator)	Z96

Microscope Lens Cleaner

A sale and effective microscope iens cleaning
solution in a dropper bottle.
30ml

Lens Paper

Non-linting, for cleaning microscope lens. 4x6 inches.

12 packs of 50 sheets.....52846001



General Lab Supplies Microscope Supplies



Bibulous Paper

Non-linting slide blotting paper. 50 sheets/book.


12 books28511007



Slide Mailer, Polypropylene

Slide mailer, polypropylene, holds 5 slides, hinged lid.

100/box.....80094376



Slide Mailer, Cardboard

Solid cardboard, holds slides for transport, with or without cover glass. Thumb groove for easy removal, cut corners identifies tip of mailer.

Holds 1 slide

2000/case.....MSL01

Holds 2 slides

1000/caseMSL02



Slide Holder Storage Box

Removable tray, holds up to 200 slides (3 inches x 1 inch).

White, each M70050W

General Lab Supplies Swabs & Applicators



Transporter Swab

Without media, rayon mini tip, individually wrapped, flexible, for urethra, ear, eye, and nasopharygeal specimens. Straight aluminum wire shaft.

50/box25800R50

Amies transport media, plastic shaft with rayon tip single swab, for use when specimens cannot be immediately plated. 50/box......4108BX

Amies Gel with Charcoal

Plastic shaft with rayon tip, single, for fastidious organisms.

50 per box.....4114BX



Tongue Blade

Smooth, splinter-free tongue blades made from Northern White Birch. Individually wrapped, sterile.

100/box.....25705



Dacron® Swabs, Sterile

Individually wrapped, wood shaft is 6 inches long.

100/box.....HD25806


Cotton Swabs, Sterile

Individually wrapped, wood shaft is 6 inches long.

100/box.....HD258061WC



General Lab Supplies Specimen Containers



SureTight[™] Specimen Sample Cups

Tight sealing lid ensures that your sample will arrive in the lab safe and dry! The patented lid is designed to securely lock in the sample. When you use SureTight[™], you can be confident that your sample is safely contained in the cup and will remain that way throughout transport to the lab! Includes tamper-evident tab, identification label, and graduation lines on sides.

Sterile, 90ml cups,

400/case......243810



Whirl-Pak® Write-On Specimen Bags

White write-on strip makes sample identification easier and more convenient. The paper-like finish virtually eliminates ink smearing and provides a surface that will accept a variety of writing instruments.

Sterile (inside), 500/box.

1oz	B01067WA
2oz	B01064WA
4oz	B01062WA
18oz	B01065WA
24oz	B01297WA
5507	R01195W/A



We have over 13,000 laboratory supply products to choose from!

General Lab Supplies Disinfectants/Cleansers, Sterilants



Alpet D2

A rapid kill, highly evaporative, no rinsing required, surface sanitizer spray approved for food contact surfaces. This isopropyl alcohol and powerful quarternary ammonium compound

eliminates 99% of detectable pathogens.

1 liter SS10003



CitriFoam® Disinfectant/ Cleaner

A powerful, phenolicbased foaming disinfectant/cleaner for use in hospitals. Broad-spectrum germicidal activity, including TB and HIV. Non-abrasive

formula won't harm surfaces. 360-degree spray nozzle works even when upside down. Mild lemon citrus scent.

12oz aerosol, 12/case.....7135



CitriGuard II® Disinfectant/ Cleaner

A new generation quaternary disinfectant that is bactericidal, virucidal, tuberculocidal, and fungicidal. EPA registered patented formula can be used in veterinary practices.

Non-acid based, biodegradable detergent safely cleans surfaces and eliminates the need for a precleaner. Ready-to-use, no dilutions needed. Contains no phenolics, glutaraldehydes, sodium hypochlorites or bleach. Active ingredient n-alkyl dimethyl benzyl ammonium chloride and ethylbenzyl ammonium chloride.

32oz spray bottle7160	
1 gallon refill bottle7142	

DisCide® Towelettes

Disinfecting towelettes impregnated with a

second-generation, dual-chain quaternary ammonium and isopropyl alcohol. The disinfectant towelette is clinically proven to kill microorganisms on hard, inanimate, nonporous surfaces, reducing the risk of cross contamination.



Destroys TB, HIV (AIDS) virus, Herpes 2 and Adenovirus 2 in one minute. Towelettes offer user-friendly non-aerosol application in the dental operatory, laboratory and surgical environments.

10.5x10.5 inches, 60/pop-up jar.....10DIS 6x7 inches, 160/pop-up jar.....60DIS

Wavicide®-01 (Sterilant)

Ready to use surface disinfectant containing 2.5% glutaraldehyde. A proven sterilant/high-level disinfectant. Requires no activator, is designed for use at room temperature. No mixing, no heating. A mild, fresh scented replacement for those pungent smelling brands, by Medical Chemical.

Wavicide® 01, sterilant, 1 gallon......01EH

General Lab Supplies Hand Soaps, Lotions, Sanitizers



AloeGuard® Antimicrobial Soap

A specially formulated antimicrobial soap with soothing, moisturizing *Aloe vera* to protect the skin's delicate balance. Provides the high antimicrobial activity necessary for

personal protection, and reduction of cross contamination. Feels and smells great while leaving your hands feeling soft and clean.

4oz flip-top bottle	.7725
18oz pump bottle	.7760
1 gallon refill	.7740

AloeSoothe[™] Moisturizing Lotion

The perfect moisturizing lotion for any healthcare professional. Because of its unique blend of lanolin, *Aloe vera* dimethicone, and other ingredients, it helps seal out harsh irritants while allowing the body's natural healing process



to occur. It also helps prevent chapping and irritation caused by frequent hand washing and use of latex gloves.

18oz pump	bottle	7750



Smart-San[™] Alcohol Mist Hand Sanitizer Spray

Ethanol based with emollients to keep hands soft, fragrance free, rated to kill 99.9%

of gram-negative and gram-positive bacteria. Requires pump, sold separately.

Six/One liter bottles	SMA0016U
Spigot (pump)	KT10011



AloeSafe[™] Hand Sanitizer Gel

Alcohol gel with *Aloe vera*, vitamin E and lavender, will not dry hands, keeps skin soft, conforms to CDC guideline for hand sanitizers. Contains no estrogenic parabens or propylene glycol (found in automobile antifreeze).

2oz pocket size	AS2HD
4 fluid oz	AS4HD
16 fluid oz	AS16HD
Wall Mounted Bottle Holder for	
16oz bottle	ASH16

Mastitis Microbiology Made Easy!

Our products are designed to aid in the rapid identification of bovine mastitis organisms and to enhance your mastitis control program.

Speed

Obtain reliable results in 24 hours for *Staphylococcus* and *Streptococcus*. Mycoplasma cultures can be detected in as little as three days.

Performance


All our mastitis media formulas are engineered for superior performance and consistency.

Availability

Hardy manufactures and sells directly to you. You can stock your entire laboratory with one purchase order. Choose from thousands of microbiology products. www.HardyDiagnostics.com

Quality

Hardy Diagnostics is known for manufacturing quality media for microbiologists since 1980.



Bovine Blood Agar Plate, Cat. no. A188



- Our Bovine Blood Agar uses bovine blood cells that are triple washed for improved sensitivity.
- Bovine Blood Agar is used for the detection of Streptococcus agalactiae and Staphylococcus aureus.
- The addition of Beta Toxin on Bovine Blood Agar aids in the identification of Streptococcus agalactiae by using a modified CAMP test.
- Staphylococcus aureus can be confirmed within one minute by using the StaphTex™ Blue Kit (Cat. no. ST50).

Bovine Blood Agar Plate

A general purpose media for the isolation and differentiation of mastitis-causing bacteria. Used to presumptively identify Staphylococcus and Streptococcus species based on colony morphology and patterns of hemolysis. 15x100mm plate, 10/pkA188

Bovine Blood Agar Plate with Esculin

Contains esculin (full strength) to differentiate group D streptococci from Streptococcus agalactiae, as S. agalactiae is not capable of esculin-hydrolysis. When esculin is hydrolyzed by organisms it forms dextrose and esculetin, which react with a compound in the media to produce a darkening or blackness around the colonies. Also known as CAMP Agar.

15x100mm plate, 10/pkA189

Bovine Blood Agar with Esculin, Modified

Bovine Blood Agar with Esculin, Modified is similar to Bovine Blood Agar with Esculin, except with half-strength esculin. Also known as washed cow blood agar.

15x100mm plate, 10/pkA186

Bovine Selective Strep Agar

For the selective isolation of *Streptococcus* species from bovine mastitis specimens. Also known as TKT Agar or modified Edwards agar. 15x100mm plate, 10/pk.....A143

Bovine Selective Staph Agar

For the selective isolation of Staphylococcus species from bovine mastitis specimens. Also known as salt agar.

15x100mm plate, 10/pk.....A157

Mastitis Biplate

Bovine Blood Agar with Esculin / MacConkey Agar 15x100mm biplate, 10/pk......J129

Mastitis Triplate

Bovine Selective Strep / MacConkey Agar / BloodAgar 15x100mm triplate, 10/pk.......J312

Mastitis Quadplate

Blood Agar / MacConkey Agar / Bovine Selective Strep Agar / Bovine Selective Staph Agar 15x100mm guadplate, 10/pk......J325



InPouch TF-Bovine

Tritrichomonas foetus in cattle is a sexually transmitted infection and is associated with aborted calves, low birth weight, or delayed calving, which can result in decreased herd vields of 8-29%. Bull infection rates are 6-8%. but have been found to be as high as 44% in some rangeland areas of the western US. Bovine Trichomoniasis is caused by the flagellated protozoan, Tritrichomonas foetus.

InPouch™ is used for the detection of nonspecific trichomonads in bovines, Trichomonas gallinae in avian, and T. suis in swine.

- Requires less hands-on time than wetmount.
- Two procedures are combined into one test (culture and identification).

10 tests	111010
100 tests	111003

Microbiology **Urinalysis**

Urine Dipstick Test

Chemstrip 10 with specific gravity (SG), urine dipstick, by Roche.

10 tests per strip, 100 strips......417145



Urine Dipstick Control Set

Normal and abnormal levels, by Quantimetrix. 4x25ml......144002

Castile Soap Towelettes

General purpose cleaning towelettes. Contain 2% coconut oil based soap solution, pH



balanced. Individual foil packets, 2 x 2.25 inches.

100 pads per box......D41900

UriStain[™]

For staining urine sediment. A stabilized modification of the Sternheimer-Malbin Urine Sediment Stain. Provides easy differentiation of red and white cells. Stains red cells pink, white cells purple, cytoplasm red-purple, yeast and fungal mycelia purple, casts pink to red, epithelial cells purple, bacteria red when alive, purple when dead, trichomonas colorless or light blue.

SureTight[®] Sample Cups

Your sample cup "dream come true!" SureTight™ represents a



new advance in container technology. The patented lid is designed to securely lock in the sample. Use SureTight™ with confidence that your sample is safely contained in the cup and will remain that way throughout transport to the lab!

- Tamper-proof sterility tab ensures sample integrity.
- Cups feature a wide 53mm opening for larger samples.
- Constructed of tough polypropylene with thick walls making them virtually indestructible.
- Ergonomic design features easy-open ridged lids.
- Paper label attached for easy sample identification. Graduation marks on the side allow for volume measurement.


Sterile, 90ml cups, 400/case.....243810



Refractometer

Handheld, comes with 3 scales for testing urine specific gravity, plasma/serum protein, and refractive index (nD), by Vee Gee. 4x4x17cm, 285g. Model CLX-1.....43027

Microbiology Urine Culture System



HardyCHROM HUrBi™ is formulated to assist in characterizing isolates by isolating gram-positive organisms on one side of the biplate and gram-negative organisms on the other side of the biplate. HUrBi™ offers chromogenic culture media to facilitate the detection and differentiation of common urinary tract microorganisms.

Logical

Total separation of organisms with grampositive on one side and gram-negative on the other side of the biplate.

Saves Money

Reduce the need for expensive automated ID cards.

Confirmatory

Identify *E. coli* and *Enterococcus* spp. with no further testing needed!

• Easy Read-Out

Distinct color reactions for each of the commonly isolated organisms. Facilitates the detection of mixed cultures.

Gram-Positive Side Gram-Negative Side

- S. aureus
- · S. epidermidis
- S. saprophyticus
- supropriy ticus
- Enterococcus spp.
- Candida spp.
- · Listeria spp.
- S. agalactiae
- Proteus, Morganella, and Providencia spp.
- · E. coli
- KES Group (Klebsiella, Enterobacter, Serratia spp.)
- spp.)
- · Citrobacter spp.
- Pseudomonas spp.

100mm biplate,

10 plates/pack.....J100


Blood Agar/MacConkey Agar Biplate

10 plates/pack......J32

Dip-Slides


CLED Agar with Polymyxin B/EMB Levine Agar. Dipslides are coated with growth media on both sides to increase surface contact with the specimen and enable the use of two different media on each slide. Use of different selective growth media permits simultaneous exposure of urine specimen to a medium selective for gram-positive species on one side, and gramnegative species on the other.

10/pack......MC500



Microbiology Skin Culture Systems (Dermatophytes)

Dermatophytes are a specialized group of parasitic fungi that most commonly include the genera *Microsporum*, *Epidermophyton* and *Trichophyton*. These fungi live only on nonviable keratin of the skin, hair and nails, and do not invade the deeper living layers of the skin, as they do not survive in these cells. They are responsible for the tinea (ringworm) infections of the skin, hair, and nails.



Derm-Duet™ is a two-section round agar plate containing DTM (Dermatophyte Test Medium) and RSM™ (Rapid Sporulation Medium™), providing the user with the simultaneous advantage of the distinctive color change of the DTM and the rapid sporulation and visual pigmentation of RSM™.

The Derm-Duet™ two-section test plate system provides for a simple comprehensive analysis of pyogenic and dermatophytic microorganisms that cause the vast majority of all infectious skin lesions seen in the clinical practice of veterinary medicine.

Derm-Duet™ II

Hardy Diagnostics now has a new Derm-Duet™ format. Individually wrapped, hermetically sealed, 12 month shelf life. 10/boxJ175


Sab-Duet[™]

DTM/Mycobiotic Agar Biplate, Mycobiotic Agar is SabDex with antibiotics for improved selection.

10/packJ850

Dermatophyte Test Medium


DTM Slant, packaged in a plastic screw-cap vial. For easy detection of dermatophytic fungi (ringworm) 20/pack.....L27





Uninoculated


Positive Reaction



DermDuet™ Individually bagged, room temperature storage, Cat. no. J175

Microbiology Staph, MRSA, Candida

Easy Identification



HardyCHROM™ Staph aureus

HardyCHROM™ Staph aureus allows for the rapid and reliable detection of *Staphylococcus aureus*. This medium contains a special chromogenic mix that allows for the differentiation of *Staphylococcus* spp.

Staphylococcus aureus can be identified as smooth, pink colored colonies on the plate. Other organisms may appear as colorless, blue or cream colonies, or will be inhibited. Staphylococcus saprophyticus will appear as turquoise colored colonies.

15x100mm plate, 10/pack......G311

HardyCHROM™ MRSA

HardyCHROM™ MRSA is a selective and differential chromogenic medium recommended for the qualitative detection of nasal colonization by methicillinresistant Staphylococcus aureus (MRSA) to aid in the prevention and control of MRSA infections in health care settings. The test is performed on anterior nares swabs from patients and healthcare workers to screen for MRSA colonization. HardyCHROM™ MRSA is not intended to diagnose MRSA infection nor to guide or monitor therapy for MRSA infections. Concomitant cultures are necessary for susceptibility testing or epidemiological typing. A negative result does not preclude MRSA nasal colonization. 15x100mm plate, 10/pack......G307



Methicillin-resistant Staphylococcus aureus (ATCC® 43300) growing on HardyCHROM MRSA, Veterinary showing magenta colonies. Incubated aerobically for 24 hours at 35° C.



C. krusei produces large rough or crenated pink to medium pink colonies.

HardyCHROM[™] **Candida**

Microbiology

Miscellaneous Culture Media

Bile Esculin Agar, BEA,

For the detection of group D streptococci including enterococci.

15x100mm plate, 10/pack......G12

Blood Agar, 5%/MacConkey Biplate

(Tryptic Soy Agar with 5% Sheep Blood / MacConkey Agar) 15x100mm biplate, 10/pack......J32

Blood Agar, 5%

For the cultivation of microorganisms. 15x100mm plate, 10/pack.A600



Brucella Agar with Hemin and Vitamin K / **LKV** (Laked Blood, Kanamycin, Vancomycin) For the cultivation of anaerobic bacteria and

the selective isolation of *Bacteroides* spp.

15x100mm biplate, 10/pack......J87

Campy CVA Agar

(Campy Agar with 5% Sheep Blood, Cefoperazone, Vancomycin, Amphotericin B) For the selective isolation of *Campylobacter jejuni* from fecal specimens.

15x100mm plate, 10/pack......A40

Chocolate Agar

(GC Agar Base with 1% hemoglobin and enrichments) For the cultivation of fastidious bacteria.

15x100mm plate, 10/pack.....E14



Columbia Blood Agar/MacConkey Biplate

(Columbia Agar with 5% Sheep Blood) For the cultivation of microorganisms, and for the isolation and differentiation of gramnegative bacteria.

15x100mm biplate,10/pack......J57

GN Broth

For the selective enrichment of Salmonella and Shigella spp. in stool specimens. 16x125mm tube, 10ml, 20/pack.....K39



EnteroScreen 4™

Single Tube Method, for Salmonella and Shigella screening from MacConkey Agar isolates from stool specimens.

16x125mmtube, 10ml, 20/pack.....L225


Microbiology Miscellaneous Culture Media

HE Agar (Hektoen Enteric) For the selective isolation of gram-negative enteric bacteria, such as Salmonella spp. and Shigella spp. 15x100mm plate, 10/pack	Rose Agar (CNA and PEA combination) (CNA and PEA Combination) For the selective isolation of gram-positive bacteria. 15x100mm plate, 10/packA66 Thioglycollate with H and K (Thioglycollate without Indicator with Hemin and Vitamin K) For the cultivation of anaerobic bacteria. 16x100mm tube, 5ml, 20/packK22
MacConkey Agar For isolation and differentiation of gramnegative bacteria. 15x100mm plate,10/pack	Saline, 0.85% For the preparation of isotonic suspensions and dilutions of microorganisms. 13x100mm tube, 2ml, 20/packR45
Mueller Hinton Agar For antimicrobial susceptibility testing of non-fastidious bacteria utilizing the disk diffusion sensitivity method. 15x150mmplate,10/packH11	SabDex C&G/DTM Biplate (Sabouraud Dextrose with Chloramphenicol and Gentamicin/Dermatophyte Test Medium) 15x100mm biplate, 10/packJ107
Mueller Hinton Agar w/Blood (Mueller Hinton Agar with 5% Sheep Blood) For disk diffusion sensitivity testing of Streptococcus spp. 15x100mm plate, 10/packA59 Mueller Hinton Agar For antimicrobial susceptibility testing of non-fastidious bacteria utilizing the disk diffusion sensitivity method.	SabDex (Sabouraud Dextrose) Agar Slant (Sabouraud Dextrose Agar) For the cultivation of fungi. 20x125mm tube, 20/packL40 Selenite Cystine Broth For the selective enrichment of Salmonella spp. Cystine enhances growth of Salmonella spp. 16x125mm, 20/packK69
15x100mm plate, 10/packG45	Thio w/o Indicator For the cultivation of microorganisms. 16x125mm tube, 10ml, 20/packK29 TSI Agar Slant, (Triple Sugar Iron Agar) For the
Mycoplasma under the microscope.	identification of gram-negative bacilli. 13x100mm tube, 20/packR32
Mycoplasma Agar For the detection of <i>Mycoplasma</i> spp. associated with bovine mastitis. 15x100mm plate, 10/packG102	Urea Agar Slant (Christensen's Urea) For the detection of urease production in bacteria. 16x100mm tube, 20/packL65

Microbiology Reagents

Bacitracin, HardyDisk™ For the differentiation of group A strep (Streptococcus pyogenes), not for sensitivity tests. Use with dispenser number 260457, by Hardy Diagnostics. 0.04 units each, 50 disks/cartridge	BIEF Water W
Catalase Reagent Hydrogen peroxide, 3%. 15ml dropper bottleZ62	
CoaguStaph™	OxiDrops [™]
Rabbit plasma with EDTA. For the detetion of coagulase enzyme in <i>Staphylococcus aureus</i> .	Oxidase reagent drops. For testing oxidase reaction.
6x20mlZ020	7ml, dropper bottleZ119
India Ink For capsule staining of <i>Cryptococcus</i> . 15ml dropper bottleZ64	OxiStrips™ Oxidase test, paper strips for testing the oxidase
Indole DMACA spot reagent. 15ml (cinnamaldehyde)Z65	reaction. 25 paper strips/pkgZ93
Nitrate Reagent A	Potassium
Sulfanilic Acid, 15mlZ71	Hydroxide
	KOH, 10% 15ml dropper bottleZ77
Nitrate Reagent B Naphthylamine,	13mi dropper bottle277
15mlZ72	PYR Test Kit
Nitrate Reagent C	A rapid test for group A strep and Enterococcus.
Zinc Dust,	25 tests, with 5ml developerZ75
15mlZ73	100 tests, with 20ml developerZ175
Optochin, HardyDisk™	Voges-Proskauer Reagent A,
Optochin, HardyDisk [™] , for the differentiation of alpha-hemolytic <i>Streptococcus pneumoniae</i> ,	Naphthol, 15ml dropper bottleZ91
use with dispenser number 260457,	
by Hardy Diagnostics. 50 disks/cartridgeZ7011	Voges-Proskauer Reagent B KOH, 40%
	15ml dropper bottleZ92

Microbiology Rapid Tests



StaphTEX[™] Blue Kit

Produces clear and reliable results when testing for *Staphylococcus aureus*. The plasma coated blue latex beads create the easiest read-out ever! Positives rapidly agglutinate the blue beads on a white background, while negatives remain consistently smooth. The StaphTEX™ kit tests for clumping factor and protein A found in *S. aureus*. Positives begin to clump immediately. The entire procedure takes less than one minute. Liquid controls, test cards, and mixing sticks are included.

50 testsST50
200 testsST200



StrepPRO™ Grouping Kit,

A latex agglutination test for the rapid identification of beta-hemolytic streps; group A, B, C, D, F, and G. All reagents are liquid, ready to use, color coded for ease of identification and simple to use. Compact container saves valuable refrigerator space. Rapid, accurate, easy to interpret positive and negative results. A complete procedure takes only 2 minutes. Kit contains reaction cards and mixing sticks. 60 tests.......PL030HD

MRSA Latex Test for PBP2'

A latex slide agglutination test for the detection of penicillin binding protein 2' (found only in MRSA) from isolates of *Staphylococcus* spp., by Denka Seiken. 50 tests......DR900A


Quality Control Microorganisms







Call (800)266-2222 for a FREE complete catalog.





Microbiology Susceptibility Testing



Antibiotic Disks: HardyDisk™ AST

Now you have a competitive choice!

Hardy Diagnostics offers an extensive selection of antibiotic sensitivity disks for the disk diffusion (Kirby-Bauer) test.

New Packaging!

HardyDisk[™] AST cartridges are packaged in an all-new canister that features a desiccant built into the lid.

Features and Benefits

- A complete line of all commonly used antibiotics.
- Spring-loaded in cartridges of 50 disks for smooth and reliable dispensing.
- · Available in single cartridges or packs of five.
- Compatible with the BD BBL™ disk dispensers (6, 8, or 12-place).
- Can also be used with the metal single disk dispenser (Cat. no. 260457).



Amikacin, AN-30 30µg, 50 disks/cartridge, 1 cartridge/canister Z8011 30µg, 50 disks/cartridge, 5 cartridges/canister.....Z8015 Amoxicillin/Clavulanic Acid, Augmentin, AmC-30 20µg/10µg, 50 disks/cartridge, 1 cartridge/canister Z8021 20µg/10µg, 50 disks/cartridge, 5 cartridges/canister..... Z8025 Ampicillin/Sulbactam, SAM-20 10μg/10μg, 50 disks/cartridge, 1 cartridge/canister......Z8041 10μg/10μg, 50 disks/cartridge, 5 cartridges/canister.....Z8045 Ampicillin, AM-10 10µg, 50 disks/cartridge, 1 cartridge/canister.....Z8031 10µg, 50 disks/cartridge, 5 cartridges/canister.....Z8035 Azithromycin, AZM-15 15µg, 50 disks/cartridge, 1 cartridge/canister.....Z8051 15µg, 50 disks/cartridge, 5 cartridges/canister.....Z8055 Bacitracin, B-10 10µg, 50 disks/cartridge, 1 cartridge/canister.......Z9271 10µg, 50 disks/cartridge, 5 cartridges/canister.....Z9275 Carbenicillin, CB-100 100µq, 50 disks/cartridge, 1 cartridge/canister.....Z8111 100µg, 50 disks/cartridge, 5 cartridges/canister.....Z8115 Cefaclor, CEC-30 30µg, 50 disks/cartridge, 1 cartridge/canister.....Z8121 30µg, 50 disks/cartridge, 5 cartridges/canister.....Z8125 Cefazolin, CZ-30 30µg, 50 disks/cartridge, 1 cartridge/canister Z8141 30µg, 50 disks/cartridge, 5 cartridges/canister.....Z8145

5µg, 50 disks/cartridge, 1 cartridge/canister......Z8181 5µg, 50 disks/cartridge, 5 cartridges/canister.....Z8185


Cefotaxime, CTX-30
30μg, 50 disks/cartridge, 1 cartridge/canisterZ8221
30μg, 50 disks/cartridge, 5 cartridges/canisterZ8225
Cefotetan, CTT-30
30μg, 50 disks/cartridge, 1 cartridge/canister Z8231
$30\mu g, 50\ disks/cartridge, 5\ cartridges/can isterZ8235$
Cefoxitin, FOX-30
30μg, 50 disks/cartridge, 1 cartridge/canisterZ8241
30μg, 50 disks/cartridge, 5 cartridges/canisterZ8245
Ceftazidime, CAZ-30
30μg, 50 disks/cartridge, 1 cartridge/canisterZ8271
30μg, 50 disks/cartridge, 5 cartridges/canisterZ8275
Ceftaroline, CPT-30
30μg, 50 disks/cartridge, 1 cartridge/canisterZ9391
30μg, 50 disks/cartridge, 5 cartridges/canisterZ9395
Ceftriaxone, CRO-30
30μg, 50 disks/cartridge, 1 cartridge/canisterZ8311
30μg, 50 disks/cartridge, 5 cartridges/canisterZ8315
Cefuroxime, CXM-30
30μg, 50 disks/cartridge, 1 cartridge/canisterZ8321
30µg, 50 disks/cartridge, 5 cartridges/canisterZ8325
Cephalothin, CF-30
30μg, 50 disks/cartridge, 1 cartridge/canisterZ8331
30µg, 50 disks/cartridge, 5 cartridges/canisterZ8335
Chloramphenicol, C-30
30µg, 50 disks/cartridge, 1 cartridge/canisterZ8341
30µg, 50 disks/cartridge, 5 cartridges/canisterZ8345
Cinoxacin, CIN-100
100µg, 50 disks/cartridge, 1 cartridge/canisterZ8351
100 µg, 50 disks/cartridge, 5 cartridges/canisterZ8355

Cefixime, CFM-5

Microbiology Susceptibility Testing

Ciprofloxacin, CIP-5	Oxacillin, OX-1
5μg, 50 disks/cartridge, 1 cartridge/canister	1μg, 50 disks/cartridge, 1 cartridge/canisterZ8851
5μg, 50 disks/cartridge, 5 cartridges/canisterZ8365	1μg, 50 disks/cartridge, 5 cartridges/canisterZ8855
Clarithromycin, CLR-15	Penicillin, P-10
15µg, 50 disks/cartridge, 1 cartridge/canisterZ8371	10μg, 50 disks/cartridge, 1 cartridge/canisterZ8881
15µg, 50 disks/cartridge, 5 cartridges/canisterZ8375	10μg, 50 disks/cartridge, 5 cartridges/canisterZ8885
Clindamycin, CC-2	Piperacillin, PIP-100
2μg, 50 disks/cartridge, 1 cartridge/canister	100μg, 50 disks/cartridge, 1 cartridge/canisterZ8901
2μg, 50 disks/cartridge, 5 cartridges/canisterZ8395	100μg, 50 disks/cartridge, 5 cartridges/canisterZ8905
Colistin Sulfate, CT	Piperacillin/Tazobactum TZP-110
10μg, 50 disks/cartridge, 1 cartridge/canisterZ8411	100/10μg, 50 disks/cartridge, 1 cartridge/canisterZ8911
Doripenem, DOR-10	100/10μg, 50 disks/cartridge, 5 cartridges/canisterZ8915
10μg, 50 disks/cartridge, 1 cartridge/canisterZ9371	Streptomycin, S-10
10μg, 50 disks/cartridge, 5 cartridges/canisterZ9375	10μg, 50 disks/cartridge, 1 cartridge/canisterZ8991
Doxycycline, D-30	10μg, 50 disks/cartridge, 5 cartridges/canisterZ8995
30μg, 50 disks/cartridge, 1 cartridge/canisterZ8441	Sulfamethoxazole/Trimethoprim, SXT
30µg, 50 disks/cartridge, 5 cartridges/canisterZ8445	1.25/23.75µg, 50 disks/cartridge, 1 cartridge/canisterZ9011
Erythromycin, E-15	1.25/23.75µg, 50 disks/cartridge, 5 cartridges/canisterZ9015
15μg, 50 disks/cartridge, 1 cartridge/canisterZ8471	Sulfisoxazole, G-0.25
15μg, 50 disks/cartridge, 5 cartridges/canisterZ8475	0.25mg, 50 disks/cartridge, 1 cartridge/canisterZ9071
Gentamicin, GM-10	0.25mg, 50 disks/cartridge, 5 cartridges/canisterZ9075
10µg, 50 disks/cartridge, 1 cartridge/canisterZ8541	Telavancin, TLV-30
10µg, 50 disks/cartridge, 5 cartridges/canisterZ8545	30µg, 50 disks/cartridge, 1 cartridge/canisterZ9381
Imipenem, IPM-10	30μg, 50 disks/cartridge, 5 cartridges/canisterZ9385
10µg, 50 disks/cartridge, 1 cartridge/canisterZ8571	Tetracycline, Te-30
10μg, 50 disks/cartridge, 5 cartridges/canisterZ8575	30μg, 50 disks/cartridge, 1 cartridge/canisterZ9121
Kanamycin, K-30	30μg, 50 disks/cartridge, 5 cartridges/canisterZ9125
30µg, 50 disks/cartridge, 1 cartridge/canisterZ8601	Ticarcillin, TIC-75
30µg, 50 disks/cartridge, 5 cartridges/canisterZ8605	75µg, 50 disks/cartridge, 1 cartridge/canisterZ9141
Levofloxacin, LVX-5	75µg, 50 disks/cartridge, 5 cartridges/canisterZ9145
5μg, 50 disks/cartridge, 1 cartridge/canisterZ8611	Tigecycline, TGC-15
5μg, 50 disks/cartridge, 5 cartridges/canisterZ8615	15µg, 50 disks/cartridge, 1 cartridge/canisterZ9361
Nalidixic Acid, NA-30	Timentin, TIM-85, (Ticarcillin with Clavulanic Acid)
30µg, 50 disks/cartridge, 1 cartridge/canisterZ8731	75µg/10µg, 50 disks/cartridge, 1 cartridge/canisterZ9161
30µg, 50 disks/cartridge, 5 cartridges/canisterZ8735	75µg/10µg, 50 disks/cartridge, 5 cartridges/canisterZ9165
Nitrofurantoin, F/M-300	Tobramycin, NN-10
300µg, 50 disks/cartridge, 1 cartridge/canisterZ8791	10μg, 50 disks/cartridge, 1 cartridge/canisterZ9171
300µg, 50 disks/cartridge, 5 cartridges/canisterZ8795	10μg, 50 disks/cartridge, 5 cartridges/canisterZ9175
Norfloxacin, NOR-10	Trimethoprim, TMP-5
10µg, 50 disks/cartridge, 1 cartridge/packZ8811	5μg, 50 disks/cartridge, 1 cartridge/canisterZ9181
10µg, 50 disks/cartridge, 5 cartridges/canisterZ8815	5μg, 50 disks/cartridge, 5 cartridges/canisterZ9185
Novobiocin, NB-30	Vancomycin, Va-30
30µg, 50 disks/cartridge, 1 cartridge/canisterZ8831	30μg, 50 disks/cartridge, 1 cartridge/canisterZ9241
30µg, 50 disks/cartridge, 5 cartridges/canisterZ8835	30μg, 50 disks/cartridge, 5 cartridges/canisterZ9245
Ofloxacin, OFX-5	
5μg, 50 disks/cartridge, 1 cartridge/canisterZ8841	
5μg, 50 disks/cartridge, 5 cartridges/canisterZ8845	


Microbiology Susceptibility Testing



Mueller	Hinton	Agar
---------	--------	------

Mueller Hinton Agar

100mm plate,
10 plates/packageG45
150mm plate,
10 plates/packageH11
Saline, 0.85 % 5ml, 16x100mm tube
,
20 tubes/packageK59



Dispensers



Multi-Dispenser	
6-place, self-tamping	.260661
8-place, self-tamping	.260660
12 place self tamping	260640



McFarland Standard

#0.5 latex equivalent, turbidity standard 16x100mm tube, 8ml, each.....ML05





Wickerham Card

For turbidity reading.
2in x 3in, each.....Z08



Sterile Swabs

Individually wrapped Dacron® swabs 100/pack......HD25806

Microbiology Loops/Inoculation



Bacti-Cinerator™ IV

Loop incinerator safely sterilizes loops and needles by infrared heat in ceramic tube: no open flame or hazardous gas. Prevents infectious spatter and/or cross contamination. Ideal for use in anaerobic chambers. Complete sterililzation takes only 5 to 7 seconds at optimum sterilizing temperature of 1500°F. Base is weighted for stability and provides convenient storage for 6 needle holders, by Sherwood.

Electric Bacti-Cinerator™ IV	004002
Ceramic Heating	
Replacement Element	.001373

Disposable SpeedStreaks™


Hardy Diagnostics' line of disposable inoculating loops are available in several styles and color coded for convenience. The 1µl and 10µl loop each come in two different styles, rigid (hard) or flexible (soft), to fit a variety of streaking preferences. Certificate of Calibration authenticates the accuracy of the loops.

1μl loop, calibrated, dark green, rigid.	
25/pack, 1000/caseHS1R	



1μl loop, calibrated, light green, flexible. 25/pack, 1000/caseHS1F
10μl loop, calibrated, dark blue, rigid. 25/pack, 1000/caseHS10R
10μl loop, calibrated, light blue, flexible. 25/pack, 1000/caseHS10F
Needle, violet, rigid,

25/pack, 1000 per case.....HSND



Reusable Calibrated Nichrome Loop


Made only from high-grade nickel-chromium alloys or pure platinum alloy to ensure repeated, accurate usage.

τμι, 26 gauge, twisted nichronie wife wi	un
aluminum handle, each	.3730
•	
10μl, 26 gauge, twisted nichrome wire v	vith
aluminum handle, each	.3740

1 ul 26 gauga twisted nichromo wire with

Microbiology Spreaders /Inoculation

Microbiology Staining Supplies





Spreaders

L- and T-shaped, sterile, disposable spreaders are an economical alternative to bending glass rods or pipets for spreading microbiological samples on agar surfaces. They are supplied sterile: no need to flame sterilize or autoclave. Unlike glass spreaders, they are non-breakable and do not pose a safety hazard. L-shaped design for convenient and "lazy" spreading. Just turn the Petri dish 360° to provide a smooth and even sample distribution.

L-Shaped Spreader, 1/pack, 500 per case	174CS01
L-Shaped Spreader, 10/pack, 500 per case	174CS10
T-Shaped Spreader, 1/pack, 500 per case	COPTS
T-Shaped Spreader, 10/pack, 500 per case	COPTS10





Advanced[™] Kit

Weak staining is a leading cause of misidentification of bacteria. Gram-positives will retain a deep blue stain, even those that are easily decolorized. Gram-negatives stain a vivid and distinct red. The combination provides a high degree of contrast for easy and reliable interpretation. Differentiate mixed cultures. When mixed on a slide, gram-positives and negatives will be easier to distinguish due to an overall brighter stain.

4x250ml bottles......GK400A



EradaStain Stain Remover A very mild


A very mild cream cleaner/ stain remover.

Effective on most histological dyes, tinctures, and general laboratory chemicals. Removes Wright's stain and most other medical stains. May be used to clean hands, clothes, plastic, glassware, and laboratory equipment, by Cambridge Diagnostics.

Microbiology Staining Supplies

Coplin Jar

For staining slides. Unit holds five single 1 inch x 3 inches slides vertically or 10 slides back-to-back. Screw cap is white linerless polypropylene, which reduces



evaporation. Rectangular base. Holds approximately 60ml. Manufactured from sodalime glass. Approximate inside dimensions: 26mm L x 26mm W x 90mm D. Jar, coplin stain, glass, screw lid.


Each......900570

Slide Staining Rack



This stainless steel staining rack is designed for use across laboratory sinks of various sizes. Triangular slide supporting rods have a tapered top, and contact the slides only in a narrow zone. Rack is simply and quickly adjustable to fit sinks up to 53.4cm (21 inch) inside dimensions. Holds slides, fits over sink, adjustable to fit most any sink.

24 inches wide......15153821



Q-Slide[™] Gram

A convenient ready-to-stain control slide that will test both your stain solutions and your staining techniques. Positive and negative stain control organisms are affixed to the circles in each Q-Slide. Up to four patient specimens (in the squares) may be processed at the same time as the QC organisms. Q-Slide Gram, a prepared slide for QC of the Gram stain contains E. coli and Staphylococcus aureus controls.

5 slides/box...... Z302

Stat Stain[™], Rapid Wright-Giemsa

For blood smears, by Hardy Diagnostics. 16oz.....SS016

Fecal Fat Stain

Acid-Fast Bacteria Stain

Acid-Fast Stain Kit, Ziehl-Neelsen, by Medical Chemical (MCC). 3x8oz.....484K

Microbiology Staining Supplies

Parasitology Diagnostic Tests



Mayer's Albumin

Used to adhere fecal specimens to glass slides in the Modified Iron Hematoxylin/ Kinyoun or Trichrome staining process, by Hardy Diagnostics.

15ml......Z69

Fecal Occult Blood QuickVue® iFOB

Detects the presence of blood in stool specimens. Immunochemical technology virtually eliminates false positive results, by Quidel.

20 tests/kit......20194



PARA-TECT™ Giardia Antigen Detection

This microwell ELISA test is an in vitro immunoassay for the qualitative determination of Giardia antigen in feces. It is a double antibody (sandwich) ELISA using an anti-giardia antibody to capture the antigen from the stool supernatant. A second anti-giardia antibody is then added which sandwiches the captured antigen. This reaction is visualized by the addition of an anti-second antibody conjugated to peroxidase and the chromogen tetramethylbenzidine (TMB). The resulting blue color development indicates the presence of Giardia antigens being bound by the anti-giardia antibodies. Results in 45 minutes with a sensitivity of 5-10 ng/ml, by Medical Chemical.

96 tests......MCCGA96

InPouch[™] TF Feline

Tritrichomonas foetus in felines is a protozoan infection of the large intestine that can lead to prolonged and intractable diarrhea. Recent studies show that T. foetus infects approximately 31% of cats. Frequently, cats are co-infected with Giardia or, because T. foetus and Giardia look similar, they are misdiagnosed with Giardia. When a cat diagnosed with Giardia fails to respond to treatment, it is recommended that the cat be tested for T. foetus.

InPouch™TF Feline can detect a single trichomonad as compared to the wetmount procedure, which is only 80% accurate. Specialized formula allows for maximum growth of trichomonads while inhibiting the growth of bacteria and *Giardia*.

- · Suitable for in-the-field testing
- No slide preparation
- · Store reagents at room temperature
- Innovative packaging

- Unique transport and test system
- Wet-mount and culture in one system
- Two procedures are combined into one test (culture and identification)



100 tests	11	1103
5 tests	11	1107

Immuno Card Stat!™ Crypto/Giardia

A rapid immunoassay for the qualitative detection of *Cryptosporidium parvum* and *Giardia lamblia* specific antigens in aqueous extracts of fecal specimens. It is intended for *in vitro* diagnostic use as an aid in the detection of suspected *Cryptosporidium* or *Giardia* infections by professional laboratories.

30 tests......750830

SAF preserved fecal smear, by Microbiologics®. 10 per box......SL1510

CONCENTRATION

ParaCon[™] Fecal Concentration Kit

Includes filter/funnels, 15ml centrifuge tubes, caps and surfactant (Triton®X), by Hardy Diagnostics.



CONTROL SLIDES

Cryptosporidium Slide

Mycobacterium and *Cryptosporidium*. Acid fast staining controls, by Microbiologics®. 10 per box.....SL4010

Intestinal Protozoa Slide

Parasitology **Transports and Stains**





TRANSPORT

Stool Transport Vial

Empty vial for stool transport. With spork in cap. Use when preservatives are not required. 30ml, 100 per case......31005



STAINS

Trichrome Stain Starter Kit

Contains Schaudinn's Fixative, Ethanol 70%, Gomori Trichrome Stain, Acid Ethanol 90%, Ethanol 95%, Carbol Xylene, and Xylene, and 10 reusable coplin jars, by Volu-Sol. 5x16oz.....VXF019



10% Formalin

Fixative recommended for preservation of fecal material used in concentration, EIA, IFA, DFA, and wet mount procedures.

15ml, 100/case......57505

Modified (Cu) PVA Fixative

Copper sulfate PVA is used for permanent staining. The copper sulfate replaces mercuric chloride.

Single Vial

Cu PVA, 100 vials/case.....280105

Two Vial Set

Cu PVA/Formalin, 50 kits/case.....89010

Para-Pak®, C&S

Culture and Sensitivity Medium (C&S) system provides a standardized procedure for the routine collection, transportation, preservation, and culture of stool specimens for bacterial enteric pathogens.

Equipment

Incubators & Thermometers





Exact-Temp[™] Thermometers

The thermometer probe stays inside the refrigerator, the digital display stays outside: no need to open the door! Temperature probe remains in plastic bottle with insulating liquid. Thermometer is traceable to National Institute of Standards and Technology (NIST) with certificate. Min/Max alarm settings with nine foot sensor cable. Choose between Celsius and Fahrenheit, ranging from -50° to 70 °C (-58 to 158 °F). Accurate up to \pm 0.6 °C for precise measurements. Resolution: 0.1 °C. Display size is 1.75 x .75 inches for easy viewing, powered by one AAA battery (included).

Certified Temperature	Temperature Range	Cat. no.
-20.0 °C	-50 °C to 70 °C (-58 °F to 158 °F)	F80301
4.0 °C	-50 °C to 70 °C (-58 °F to 158 °F)	R80301
37.0 °C	-50 °C to 70 °C (-58 °F to 158 °F)	180301

Incubator

This compact incubator fits easily on most counter tops. The incubator holds a constant temperature to plus or minus one degree. It is made of powder coated steel and has an acrylic see-through door. Temperature range from ambient +5 °C up to 60 °C. Capacity is 0.8 cubic feet. The inside dimensions of the chamber are:

Overall dimensions: 13"x14"x16-1/4", Weight: 19 lbs.....132000



Hardy Diagnostics ships from 7 U.S. Distribution Centers:

Santa Maria, California Centralia, Washington Salt Lake City, Utah Phoenix, Arizona Dallas, Texas Springboro, Ohio Lake City, Florida ...which means you get your order





HARPIN A Culture of Service"
1430 West McCoy Lane
DIAGNOSTICS Santa Maria, CA 93455 1430 West McCoy Lane Santa Maria, CA 93455